


NOTES. 1. DEBURR OR BREAK ALL EDGES.
2. ▽ INDICATES FEATURES TO BE
CHECKED AT INCOMING Q.C..



Both .38 dms must be
maintained .53 deep
both ends of part.
Per R. Murray-BCT
DUE 4/26/01

CAD GENERATED DWG
850845-000.prt

REV	REVISION DESCRIPTION	BY	DATE	ECO *
A	ORIGINAL RELEASE	MA	12/15/98	58374

Dated Print 4/26/01

THESE DIMENSIONS AND SPECIFICATIONS ARE THE PROPERTY OF COBE BCT, INC. AND SHALL NOT BE REPRODUCED, COPIED, OR USED FOR MANUFACTURE OR SALE OF ANY PART WITHOUT THE WRITTEN PERMISSION OF COBE BCT, INC.		COBE BCT, Inc. Blood Component Technology 1201 Oak Street Lakewood, Colorado 80215-4407 USA © Gemini Company	
TOLERANCES AND SYMBOLS PER ANSI Y14.5M-1992 UNLESS OTHERWISE SPECIFIED DIMS. FINISH		TITLE	
.XX .00 ANGLES .XX1 .02 ± 2 .XXX .010 SURF. FINISH .XXXX .0000 125		MATERIAL STAINLESS STEEL	
DATE 12/15/98		DATE 12/15/98	
ORIGINATOR M2		DATE 12/15/98	
SIZE B		DATE 12/15/98	
SCALE 1:1		DATE 12/15/98	
DWG 850845-000.prt		DATE 12/15/98	
REV A		DATE 12/15/98	
ECO 11		DATE 12/15/98	

TOOL: JAW ADJUSTING

850845-000

11/1

408

1/24/94